

Prototype: CNN + Grad-CAM on CBIS-DDSM subset

Import Libraries

```
# %%  
# Prototype: CNN + Grad-CAM on CBIS-DDSM subset  
  
import os  
import sys  
import numpy as np  
import pandas as pd  
import matplotlib.pyplot as plt  
  
import torch  
import torch.nn as nn  
from torch.utils.data import DataLoader  
from torchvision import models  
  
# Make project root importable  
ROOT_DIR = os.path.abspath(os.path.join(".."))  
if ROOT_DIR not in sys.path:  
    sys.path.append(ROOT_DIR)  
  
from utils.dataset import CBISDicomDataset # you'll create/adjust this  
from utils.transforms import get_preprocessing_transforms  
from utils.gradcam import GradCAM  
  
device = "cuda" if torch.cuda.is_available() else "cpu"  
device  
  
'cpu'
```

Paths

```
RAW_ROOT = os.path.join(ROOT_DIR, "data", "raw", "CBIS-DDSM-All-doiJNLP-zzWs5zfZ", "cbis_ddsm_raw")  
METADATA_MASS_TRAIN = os.path.join(ROOT_DIR, "data", "raw", "mass_case_description_train_set.csv")  
METADATA_MASS_TEST = os.path.join(ROOT_DIR, "data", "raw", "mass_case_description_test_set.csv")  
  
print(RAW_ROOT)  
print(METADATA_MASS_TRAIN)  
  
/Users/gibionmakiwa/MyDocuments/610MtSundanceBayW/University_of_London_Computer_Science_BSc/  
/Users/gibionmakiwa/MyDocuments/610MtSundanceBayW/University_of_London_Computer_Science_BSc/
```

Load metadata and build a small prototype subset (e.g., 200 cases)

```
CSV_PATH = os.path.join(
    ROOT_DIR,
    "data",
    "raw",
    "CBIS-DDSM-All-doiJNLP-zzWs5zfZ",
    "metadata",
    "mass_train_with_dicom_paths.csv"
)

meta = pd.read_csv(CSV_PATH)

# Keep only rows that have a valid DICOM path
meta = meta[meta["actual_dicom_path"].notna()].reset_index(drop=True)

# Sample 200 for prototype
meta_subset = meta.sample(n=200, random_state=42).reset_index(drop=True)

# Ensure label is binary
meta_subset["label"] = (meta_subset["pathology"] == "MALIGNANT").astype(int)

meta_subset.head()
```

	patient_id	breast_density	left or right	breast image view	abnormality id	\
0	P_00976	1	LEFT	MLO	1	
1	P_01394	2	RIGHT	MLO	1	
2	P_00871	2	LEFT	CC	1	
3	P_00068	3	RIGHT	CC	1	
4	P_01849	1	RIGHT	MLO	1	

	abnormality type	mass shape	mass margins	assessment	\
0	mass	IRREGULAR	OBSCURED	0	
1	mass	IRREGULAR	SPICULATED	5	
2	mass	IRREGULAR	ILL_DEFINED	4	
3	mass	IRREGULAR	SPICULATED	4	
4	mass	LYMPH_NODE	NaN	2	

	pathology	subtlety	\
0	BENIGN	2	
1	MALIGNANT	5	
2	MALIGNANT	3	
3	MALIGNANT	4	
4	BENIGN_WITHOUT_CALLBACK	3	

```

                                image file path \
0 Mass-Training_P_00976_LEFT_MLO/1.3.6.1.4.1.959...
1 Mass-Training_P_01394_RIGHT_MLO/1.3.6.1.4.1.95...
2 Mass-Training_P_00871_LEFT_CC/1.3.6.1.4.1.9590...
3 Mass-Training_P_00068_RIGHT_CC/1.3.6.1.4.1.959...
4 Mass-Training_P_01849_RIGHT_MLO/1.3.6.1.4.1.95...

                                cropped image file path \
0 Mass-Training_P_00976_LEFT_MLO_1/1.3.6.1.4.1.9...
1 Mass-Training_P_01394_RIGHT_MLO_1/1.3.6.1.4.1....
2 Mass-Training_P_00871_LEFT_CC_1/1.3.6.1.4.1.95...
3 Mass-Training_P_00068_RIGHT_CC_1/1.3.6.1.4.1.9...
4 Mass-Training_P_01849_RIGHT_MLO_1/1.3.6.1.4.1....

                                ROI mask file path \
0 Mass-Training_P_00976_LEFT_MLO_1/1.3.6.1.4.1.9...
1 Mass-Training_P_01394_RIGHT_MLO_1/1.3.6.1.4.1....
2 Mass-Training_P_00871_LEFT_CC_1/1.3.6.1.4.1.95...
3 Mass-Training_P_00068_RIGHT_CC_1/1.3.6.1.4.1.9...
4 Mass-Training_P_01849_RIGHT_MLO_1/1.3.6.1.4.1....

                                tcia_folder \
0 Mass-Training_P_00976_LEFT_MLO_1
1 Mass-Training_P_01394_RIGHT_MLO_1
2 Mass-Training_P_00871_LEFT_CC
3 Mass-Training_P_00068_RIGHT_CC_1
4 Mass-Training_P_01849_RIGHT_MLO

                                actual_dicom_path  label
0 /Users/gibionmakiwa/MyDocuments/610MtSundanceB...      0
1 /Users/gibionmakiwa/MyDocuments/610MtSundanceB...      1
2 /Users/gibionmakiwa/MyDocuments/610MtSundanceB...      1
3 /Users/gibionmakiwa/MyDocuments/610MtSundanceB...      1
4 /Users/gibionmakiwa/MyDocuments/610MtSundanceB...      0

```

Inspect how DICOM paths are encoded in metadata

```

# %%
# Inspect how DICOM paths are encoded in metadata

meta_subset[["patient_id", "image file path", "label"]].head()

  patient_id                                image file path  label
0    P_00976  Mass-Training_P_00976_LEFT_MLO/1.3.6.1.4.1.959...      0
1    P_01394  Mass-Training_P_01394_RIGHT_MLO/1.3.6.1.4.1.95...      1
2    P_00871  Mass-Training_P_00871_LEFT_CC/1.3.6.1.4.1.9590...      1

```

3	P_00068	Mass-Training_P_00068_RIGHT_CC/1.3.6.1.4.1.959...	1
4	P_01849	Mass-Training_P_01849_RIGHT_MLO/1.3.6.1.4.1.95...	0

Create PyTorch Dataset using utils/dataset.py

```
# %%
# Create PyTorch Dataset using your utils/dataset.py

from torch.utils.data import random_split

transform = get_preprocessing_transforms() # from utils/transforms.py

full_dataset = CBISDicomDataset(
    metadata=meta_subset,
    load_masks=True
)

len(full_dataset)

<utils.dataset.CBISDicomDataset object at 0x12936acd0>
```

Split into train/val for the prototype

```
# %%
# Split into train/val for the prototype

train_size = int(0.8 * len(full_dataset))
val_size = len(full_dataset) - train_size
train_dataset, val_dataset = random_split(full_dataset, [train_size, val_size])

train_loader = DataLoader(train_dataset, batch_size=16, shuffle=True, num_workers=2)
val_loader = DataLoader(val_dataset, batch_size=16, shuffle=False, num_workers=2)

len(train_dataset), len(val_dataset)

(160, 40)
```

Visual sanity check: show one sample

```
# %%
# Visual sanity check: show one sample

img, mask, label = val_dataset[0]

input_tensor = img.unsqueeze(0).to(device)
cam = gradcam(input_tensor)
```

```

img_np = img.permute(1,2,0).numpy()
mask_np = mask.squeeze().numpy()

plt.figure(figsize=(12,4))

plt.subplot(1,3,1)
plt.imshow(img_np, cmap="gray")
plt.title("Original")

plt.subplot(1,3,2)
plt.imshow(mask_np, cmap="gray")
plt.title("ROI Mask")

plt.subplot(1,3,3)
plt.imshow(img_np, cmap="gray")
plt.imshow(cam, cmap="jet", alpha=0.4)
plt.title("Grad-CAM Overlay")

plt.show()

```

```

error                                     Traceback (most recent call last)

```

```

Cell In[7], line 4

```

```

     1 # %%
     2 # Visual sanity check: show one sample
---->  4 img, label = full_dataset[0]
       5 img_np = np.transpose(img.numpy(), (1, 2, 0))
       6 plt.imshow(img_np, cmap="gray")

```

```

File ~/MyDocuments/610MtSundanceBayW/University_of_London_Computer_Science_BSc/CM3070 Comput

```

```

    26 img = np.stack([img, img, img], axis=-1)
    28 if self.transform:
---->  29     img = self.transform(img)
    31 return img, label

```

```

File /usr/local/lib/python3.8/site-packages/torchvision/transforms/transforms.py:95, in Comp

```

```

    93 def __call__(self, img):
    94     for t in self.transforms:
---->  95         img = t(img)
    96     return img

```

```

File ~/MyDocuments/610MtSundanceBayW/University_of_London_Computer_Science_BSc/CM3070 Comput

```

```

     8 if img.ndim == 2:
     9     img = cv2.cvtColor(img, cv2.COLOR_GRAY2RGB)
---->  11 img = cv2.equalizeHist(img[:, :, 0])
    12 img = cv2.cvtColor(img, cv2.COLOR_GRAY2RGB)

```

```

13 img = cv2.resize(img, (224, 224))

error: OpenCV(4.12.0) /Users/runner/work/opencv-python/opencv-python/opencv/modules/imgproc/

```

Build and configure ResNet50

```

# %%
# Build and configure ResNet50

resnet = models.resnet50(weights=models.ResNet50_Weights.IMAGENET1K_V2)
num_features = resnet.fc.in_features
resnet.fc = nn.Linear(num_features, 1)

resnet = resnet.to(device)

criterion = nn.BCEWithLogitsLoss()
optimizer = torch.optim.Adam(resnet.parameters(), lr=1e-4)

```

Training loop (prototype: ~5–10 epochs)

```

# %%
# Training loop (prototype: ~5-10 epochs)

from tqdm.auto import tqdm

def train_one_epoch(model, loader, optimizer, criterion, device):
    model.train()
    running_loss = 0.0

    for imgs, labels in tqdm(loader, desc="Train", leave=False):
        imgs = imgs.to(device)
        labels = labels.float().to(device)

        optimizer.zero_grad()
        outputs = model(imgs).squeeze(1)
        loss = criterion(outputs, labels)
        loss.backward()
        optimizer.step()

        running_loss += loss.item() * imgs.size(0)

    return running_loss / len(loader.dataset)

def evaluate(model, loader, criterion, device):

```

```

model.eval()
running_loss = 0.0
all_outputs = []
all_labels = []

with torch.no_grad():
    for imgs, labels in tqdm(loader, desc="Val", leave=False):
        imgs = imgs.to(device)
        labels = labels.float().to(device)

        outputs = model(imgs).squeeze(1)
        loss = criterion(outputs, labels)

        running_loss += loss.item() * imgs.size(0)
        all_outputs.append(outputs.cpu())
        all_labels.append(labels.cpu())

all_outputs = torch.cat(all_outputs)
all_labels = torch.cat(all_labels)

probs = torch.sigmoid(all_outputs)
preds = (probs > 0.5).int()
acc = (preds == all_labels.int()).float().mean().item()

return running_loss / len(loader.dataset), acc

EPOCHS = 10
for epoch in range(EPOCHS):
    train_loss = train_one_epoch(resnet, train_loader, optimizer, criterion, device)
    val_loss, val_acc = evaluate(resnet, val_loader, criterion, device)
    print(f"Epoch {epoch+1}/{EPOCHS} | Train Loss: {train_loss:.4f} | Val Loss: {val_loss:.4f} | Val Acc: {val_acc:.4f}")

```

Save prototype model

```

## # %%
# Save prototype model

os.makedirs(os.path.join(ROOT_DIR, "models"), exist_ok=True)
MODEL_PATH = os.path.join(ROOT_DIR, "models", "resnet50_prototype.pth")
torch.save(resnet.state_dict(), MODEL_PATH)
MODEL_PATH

```

Grad-CAM: use utils/gradcam.py

```
# %%  
# Grad-CAM: use your utils/gradcam.py  
  
# We assume utils/gradcam.py defines a GradCAM class similar to:  
# GradCAM(model, target_layer)  
  
target_layer = resnet.layer4[-1]  
gradcam = GradCAM(resnet, target_layer, device=device)
```

Pick a validation sample and generate Grad-CAM

```
# %%  
# Pick a validation sample and generate Grad-CAM  
  
resnet.eval()  
  
img, label = val_dataset[0]  
input_tensor = img.unsqueeze(0).to(device)  
  
cam = gradcam(input_tensor) # expected shape: (H, W) normalised 0-1  
  
# Convert original image back to HxW or HxWx3 for overlay  
img_np = np.transpose(img.numpy(), (1, 2, 0))  
  
heatmap = plt.cm.jet(cam)[..., :3] # RGB  
overlay = 0.4 * heatmap + 0.6 * img_np  
  
plt.figure(figsize=(12, 4))  
plt.subplot(1, 3, 1)  
plt.imshow(img_np, cmap="gray")  
plt.title(f"Original (label={int(label)})")  
plt.axis("off")  
  
plt.subplot(1, 3, 2)  
plt.imshow(cam, cmap="jet")  
plt.title("Grad-CAM")  
plt.axis("off")  
  
plt.subplot(1, 3, 3)  
plt.imshow(overlay)  
plt.title("Overlay")  
plt.axis("off")  
  
plt.tight_layout()
```



```
plt.show()
```

Loop over several samples and save Grad-CAM overlays to results/heatmaps

```
# %%
```

```
# (Optional) Loop over several samples and save Grad-CAM overlays to results/heatmaps
```

```
HEATMAP_DIR = os.path.join(ROOT_DIR, "results", "heatmaps")  
os.makedirs(HEATMAP_DIR, exist_ok=True)
```

```
for i in range(10):
```

```
    img, label = val_dataset[i]
```

```
    input_tensor = img.unsqueeze(0).to(device)
```

```
    cam = gradcam(input_tensor)
```

```
    img_np = np.transpose(img.numpy(), (1, 2, 0))
```

```
    heatmap = plt.cm.jet(cam)[..., :3]
```

```
    overlay = 0.4 * heatmap + 0.6 * img_np
```

```
    plt.imsave(os.path.join(HEATMAP_DIR, f"sample_{i}_overlay_label_{int(label)}.png"), overlay)
```